TH WORLD STERILIZATION CONGRESS 2018

XIII INTERNATIONAL STERILIZATION CONGRESS AND HOSPITAL DISINFECTION

OCTOBER 31 TO NOVEMBER 3 2018 WORLD TRADE CENTER MEXICO CITY

SCIENTIFIC PROGRAM









Automation in the CSSD Benefits and opportunities


Consultent Gitte Antonisen Antonisen Consult



The Danish expertgroup recommendations

- Centralizing of reprocessing
- Case cart deliveries
- Use of automation
- Standardization of instrument trays
- Increased use of disposable instruments







Pilot projects were automation and use of robotics were in focus.

- The Defu Stepp project was running over three years from 2010.
 15 companies joined as partners.
- The goal was the development of a fully automatic CSSD. The project had four focus area:


Sub-projects:


- Overall production planning in the CSSD
- Automation of sub-processes in the unclean area in the CSSD
- Development of autoclavable case cart trolley which also form the sterile barrier
- Automatic attachment and removal of instrument containers (cabinet washing machine)


1.7 million new industrial robots by 2020

IFR
International
Federation of
Reportics


Double-digit average annual increase












Automation CSSD - opportunities

SUITABLE FOR AUTOMATION	TECHNOLOGIES
Transportation of goods	Transportation systems AGVs and transport cabinets (shuttler) Storage cranes and lifts Robots
Storage facilities	High-bay racking system, i.e. automatic mini load storage Storage robot
Lift, push, stabling, loading/unloading, open/close, handling etc.	Robots AGVs (Automated Guided Vehicles)
Control systems (recognizability and weight)	Vision- and weight systems
Control of flow, functions and locations	MES (Manufacturing Execution Systems) WMS (Warehouse Management Systems)

2 Centralized, automated sterilization centers



Conditions and requirements

- The two central sterile units must be identical
- Design and work flow must be uniform
- Workflow to automate as far as possible and must be identical for the two sterilization centers
- Maintenance of equipment should be coordinated as much as possible
- · Use of standard instruments

1. AUTOMATIC TRANSPORT OF GOODS


The "good old days"...




Now or future

- AGVs = Automated Guided Vehicles
- Roller conveyer transfer cabinets (shuttles)
- · Transport systems
- Transport trolleys
- · Autoclave racks
- Washing racks with instrument containers etc.
- · Washing stands with instruments
- Instrument containers
- Transport containers
- Instrument net






AGV – Boxes, instrument containers, racks – and placing item on/off washing stands



AGV – Boxes, instrument containers, racks – and placing item on/off washing stands




Roller conveyer - transfer cabinets (shuttles)



2. STORAGE FACILITIES

The "good old days"...






3. HANDLING, LOADING & ONLOADING

The "good old days"...



Now or future




3. HANDLING, LOADING & ONLOADING

The "good old days"...





4. CONTROL SYSTEM & FLOWS



CHALLENGES

- Hygienic design
- Space for cleaning




CHALLENGES

- Hygienic design
- Space for cleaning




CHALLENGES


- Hygienic design
- Space for cleaning





The Intelligent reprocessing unit-Automated Packaging and Inspection (DISAPI)






Fully automatic cleaning control

- Based on multi-spectral vision, and automation technology
- Identify debris on instruments
- Evaluation of the instrument is clean, unclean or in need of further control
- Potential for detection of corrosion damage






Vision - recognition of instruments


 The solution uses a combination of camera and laser scanning

 The position of the instruments in the trays is dented - Robot gripper

 Gripper, which is so flexible that it can handle all the selected test instruments



Fully automatic packing table



Benefits and opportunities

KEY WORDS

- Minimization of errors increased patient safety
- Precision Documentation and Traceability
- Better working environment Ergonomics
- Increased efficiency